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THE INFORMATION REQUIRED TO ENABLE THE USER TO IDENTIFY THE PRODUCT AND ITS CONTENTS, THE DESCRIPTION OF THE PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT TECHNOLOGY, APPLIED TO

Cages / spacers are implants used in several countries, and their use, application and characteristics are enshrined in vertebral surgery. The cages are small, rectangular, hollow, usually made of titanium. These spacers are often used to restore the natural high of the intervertebral space due to a compromised disk and / or to relieve pressure on the root. Cages should have an internal space to be filled with bone graft between the two vertebrae in order to facilitate the fusion between hem (arthrodesis), which increases the stability of the system in the long term. The implants may be filled with any bone graft, and the choice is the surgeon's responsibility. The bone graft begins to grow through the existing holes in the cages walls fo solid bone (fusion) which shall hold the vertebrae together. This process is known as intervertebral fusion.

The Non-expansive Cage BM, is an anterior implantable medical product, when using the cervical cages, or posterior, when sing the lumbar spacers. The non-expansive Cages BM are composed of Cervical Cages and Lumbar Spacers The finishing of implants is made by mechanical polishing, then the product goes through a pickling process, which consists of an acid bath that removes the impurities of the manufacturing process, such as oils, greases and other kinds of manufacturing materials. Then, the product goes through a process of ionization, that is, a superficial treatment which seeks to strengthen the superficial layer of products made of titanium. The Non-Expansive Cages BM are wrapped up in properly identified plastic

ackaging. The graphical information presented below are illustrative only, and refer to the non-expansive Cages - BM

The non-expansive Cages - BM are made of Titanium Alloy according to the specifications of the standard ASTM F136. Titanium is known for its extreme chemical passivity, and excellent biocompatibility, consequently it has physical properties suitable for a good biomechanical behavior in the long term. The implants made of titanium have about 45% less density than those containing iron and cobalt in their compositions, an important factor related to the patient's comfort. Its low modulus of elasticity is onsidered to be another advantage, because it minimizes the protection against pressure, and it is transferred to the bone. This elative importance against pressure increases as the size of the implant increases. The titanium is extremely insoluble, and serves as an inert material that does not interact with the body (Martinez 2004). According to Vuillemin et al. there is a physiological saturation of titanium in the body, therefore the possibility of interaction with any additional soluble titanium is unlikely, and allergic or toxic reactions were not described. The modulus of elasticity of titanium is relatively low (E = 110 000N/mm2). It is the lowest of the metals. Its elasticity is twice the cobalt-chromium alloy and four times less elastic than the ortical bone. It is a material of extreme hardness and this can only be shared with the ceramics

he non-expansive Cages - BM will be distributed on a uniform basis, according to the following table: Product Name / Description Size (width x length x thickness)

8450-05-012	Cervical Cage, 0 degree - Titanium	14mm x 12mm x 5mm
8450-06-012	Cervical Cage, 0 degree - Titanium	14mm x 12mm x 6mm
8450-07-012	Cervical Cage, 0 degree - Titanium	14mm x 12mm x 7mm
8450-08-012	Cervical Cage, 0 degree - Titanium	14mm x 12mm x 8mm
8450-05-014	Cervical Cage, 0 degree - Titanium	16mm x 14mm x 5mm
8450-06-014	Cervical Cage, 0 degree - Titanium	16mm x 14mm x 6mm
8450-07-014	Cervical Cage, 0 degree - Titanium	16mm x 14mm x 7mm
8450-08-014	Cervical Cage, 0 degree - Titanium	16mm x 14mm x 8mm
8451-05-012	Cervical Cage, 4 degree - Titanium	14mm x 12mm x 5mm
8451-06-012	Cervical Cage, 4 degree - Titanium	14mm x 12mm x 6mm
8451-07-012	Cervical Cage, 4 degree - Titanium	14mm x 12mm x 7mm
8451-08-012	Cervical Cage, 4 degree - Titanium	14mm x 12mm x 8mm
8451-05-014	Cervical Cage, 4 degree - Titanium	16mm x 14mm x 5mm
8451-06-014	Cervical Cage, 4 degree - Titanium	16mm x 14mm x 6mm
8451-07-014	Cervical Cage, 4 degree - Titanium	16mm x 14mm x 7mm
8451-08-014	Cervical Cage, 4 degree - Titanium	16mm x 14mm x 8mm
8900-03-000	Lumbar Spacer, 0 degree - Titanium	11mm x 25mm x 6mm
8900-07-000	Lumbar Spacer, 0 degree - Titanium	11mm x 25mm x 7mm
8900-08-000	Lumbar Spacer, 0 degree - Titanium	11mm x 25mm x 8mm
8900-09-000	Lumbar Spacer, 0 degree - Titanium	11mm x 25mm x 9mm

8900-10-000 Lumbar Spacer, 0 degree - Titanium 8900-11-000 Lumbar Spacer, 0 degree – Titanium 8900-12-000 Lumbar Spacer, 0 degree – Titanium Lumbar Spacer, 0 degree – Titanium Lumbar Spacer, 0 degree – Titanium 8900-06-001 Lumbar Spacer, 4 degree – Titanium 11mm x 25mm x 6mm Lumbar Spacer, 4 degree – Titanium Lumbar Spacer, 4 degree – Titanium Lumbar Spacer, 4 degree – Titanium 8900-09-001 Lumbar Spacer, 4 degree – Titanium 11mm x 25mm x 9mm 11mm x 25mm x 10mm 8900-13-001 Lumbar Spacer, 4 degree – Titanium 8900-14-001 Lumbar Spacer, 4 degree - Titanium 11mm x 25mm x 14mm

8585-00-000

A	ILENTION	
In	order to use the	non-expansive Cages it is necessary the use of the below specified instrumer
Th	ne instruments r	nust be purchased separately from Cages, because they are not part of this pro
1	8640-00-000	Rasp
1	8641-00-000	Graft Cutter/Trimmer
1	8642-00-000	Curved Impactor (Lateral) 1 8643-00-000 – Curved Impactor (Superior)
1	8644-70-000	Parallel Distractor
1	8644-90-000	Parallel Distractor EX
1	8645-00-000	Universal Modeler
1	8646-00-000	Positioner wrench(inseter) for small implants
1	8647-00-000	Positioner wrench (inserter) for large implants

8648-00-000 Positioner of small implants 8649-00-000 small Impactor Long Impactor Cutting Forceps Screw Tweezer 6157-00-000

Positioner/Impactor Key

Quick-release "T" handle Quick-release curette Bone Graft Impactor 8583-05-000 Graft Impactor Base 8587-00-000 Pin-fixing Key 8582-00-000 8582-05-000 Trephine Internal Stem

8581-00-000 8580-00-000 Puricture 8580-00-000 Cervical Sc 8584-00-000 Distractor 8588-05-012 Cervical Cage 0 º 12 x 5mm (Test) 8588-06-012 Cervical Cage 0 ° 12 x 6mm (Test)
8588-05-014 Cervical Cage 0 ° 12 x 7mm (Test)
8588-05-014 Cervical Cage 0 ° 14 x 5mm (Test)

8588-06-014 Cervical Cage 0° 14 x 6mm (Test 8588-07-014 Cervical Cage 0° 14 x 7mm (Test **8589-05-012** Cervical Cage 4º 12 x 5mm (Test) **8589-06-012** Cervical Cage 4º 12 x 6mm (Test) 8589-07-012 Cervical Cage 4º 12 x 7mm (Test 8589-05-014 Cervical Cage 4° 14 x 5mm (Test)
8589-05-014 Cervical Cage 4° 14 x 6mm (Test)
8589-07-014 Cervical Cage 4° 14 x 7mm (Test)
Cervical Cage 4° 14 x 7mm (Test)

8686-00-000 Diapason for Lumbar spacers 8682-07-000 Intervertebral Spacer Ø7.0 Intervertebral Spacer Ø8. 0 8682-10-000 Intervertebral Spacer -10, 0 8682-11-000 Intervertebral Spacer Ø11.0

Nerve Tissue Spacer 18680-00-000 – Beater for lumbar spacer Positioner Impactor Spacer stem with thread Graf Impactor Key 8684-00-000 Lumbar Positioner Impactor 8685-07-000 Curette for Lumbar Spacer - 07 3683-07-000 Intervertebral Lumbar Scraper Ø7, 0 8683-08-000 Intervertebral Lumbar Scraper Ø8. 0 8683-09-000 Intervertebral Lumbar Scraper Ø9, 0 8683-10-000 Intervertebral Lumbar Scraper -10, 0

PRODUCT PRESENTATION FOR CONSUMPTION

8684-01-000

Each packaging of the implantable component has two labels containing the data of the implant used. One of the labels should be affixed to the patient's medical records, and the other is for the hospital control. The hospital must register this information in their files and give it to the patient. In the medical records, the following information is essential:

name of the implant used / surgery date / product code / batch number / name of the patient who received the implant / name of

The implants are packed in plastic envelopes containing a non-sterile label, an external label containing the implant data, an internal instruction for use, and two internal labels containing the data of implant used. The product has the following information a) Commercial Product Name b) Technical Name c) Batch Number d) Registration Number at the Ministry of Health; e) Product Code f) Quantity g) Description of the product that comes in the package h) Date of manufacture; i) Expiration Date: Indefinite j) Description of the raw material used to manufacture the product; k) Manufacturer's address/written in the plastic envelope); l) Name of Responsible technician; m) and the statement "Product not sterile."; N) Security Symbols printed in the plastic envelopes as shown below: o) "Product for Single-Use, do not reuse it." P) Before using, see instructions for use.

SPECIAL CONDITIONS FOR STORAGE, CONSERVATION AND / OR MANIPULATION OF THE PRODUCT. NECESSARY CARE WITH HANDLING AND TRANSPORT OF THE MEDICAL PRODUCT

-Do not use the product if it is damaged.
-The implants should be handled with care to avoid harm that could undermine the quality of equipment and patient safety. -Store it in a dry, airy place away from direct sunlight and bad weather.
-When transporting it, you should avoid collisions, and inadequate stacks

TESTS MADE (LOAD RESTRICTIONS / MECHANICAL RESISTANCE)

Indication, purpose or use of the medical product as indicated by the manufacturer.

The NON-EXPANSIVE CAGES - BM have the purpose or indication of use to promote a stable biomechanical structure to the spine in order to facilitate arthrodesis (vertebral fusion). The NON-EXPANSIVE CAGES - BM are not reusable.

NOTE: The implants shown here were designed specifically for the use as described above. Any other use is contraindicated.

The Cervical Cage is indicated in cases of anterior intersomatic arthrodesis. It can be used alone or with the Implant System for Spine Fixation - BM Plates and / or the Pedicular Fixation System. Its ends, top and bottom, are formed by a serrated surface

ISO 9001 ISO 13485

allows the use of the instrumen

that is designed to better attachment to the vertebral segment avoiding migration. Its interior allows the filling with bone graft. The Lumbar Spacer is indicated for lumbar intersomatic arthrodesis by posterior approach, and in case of degenerative disk disease. This implantable metal component can be used alone or with the implant system for Spine Fixation - BM plates. The Lumbar Spacer has a rectangular shape with an internal hole. Its ends, top and bottom, are formed by a serrated surface that is designed to better attachment to the lumbar segment avoiding migration. Its top and bottom parts are composed by an oblong

tear and three more holes on the smooth sides where bone graft is placed. There is a threaded orifice on its posterior part that

The Surgical techniques depend on the surgeon, which is in charge of choosing the proper method, type and size of the products being used.

The correct selection of the implants is essential for the surgery success. Scratches or cuts must be avoided when handling the implants, since these defects are hubs of tension and can be areas of crack nucleation and reduce the resistance to corre

e below the maximum values found in the compression/shear, compression and fatigue of non-expansive – BM Cages tests

Prod	Test		Average
	Compression	Compression Maximum load	25610,20N
l Cage	Compression / Shear	Compression Maximum load Shear	9830,12N
Cervical Cage	Fatigue	Load applied during 5 million cycles under a 10Hz frequency	17927,14N
	Roughness	Roughness (Ra)	0,36µm
	Compression	Compression Maximum load	24668,80N
Spacer	Compression / Shear	Compression Maximum load Shear	10310,20N
Lumbar S	Fatigue	Load applied during 5 million cycles under a 10Hz frequency	18501,46N
Ē	Roughness	Roughness (Ra)	0,472µm

ALL WARNINGS AND/OR PRECAUTIONS THAT MUST BE ADOPTED

Single-Use Product, They must be destroyed after explantation. Do not reuse the product, An orthopedic implant can be used only in one patient, only once. Even though it may seem undamaged, the success of the implant might be reduced because of ections caused by prior tensions. Selection of improper implant can cause unusual tensions and may result in

monitor the restoration of the patient through the X-ray, according to international safety standards. Note: The use of magnetic resonance imaging in patients who already have implants can damage the diagnosis. It can also cause the displacement of the implant.

Obesity has an effect on the spinal system load, and can compromise the spinal implant fixation. Excessive physical activity and

remembered that in spite of having biocompatible material, the product, when deployed, is subject to constant environmenta changes promoting the concentration of salts, acids and alkaloids and they can cause corrosion. The clinical results and

The use of this implant is contraindicated in the following situations

including: Compression test based on ASTM F2077: 2003 (Test Methods For Intervertebral Body Fusion Devices).
 Fatigue test based on ASTM F2077: 2003 (Test Methods For Intervertebral Body Fusion Devices).

subsequent fracture.

Due to its form of presentation and its characteristics, this product does not have side effects, however, it is suggested to international safety standards. Note: The use of magnetic

traumatism directed linked to the thoracic and cervical segments may jeopardize the success of the implant prematurely. We recommend that the activities are measured avoiding excessive "stress" in the operated area. The Non-expansive Cages - BM may not be under loads or levels of activity similar to normal and healthy bones.

The mixture of different biomedical materials not compatible and/or similar may raise the risk of corrosion. It must be

durability of the implant depend on the existence of a precise surgical technique. Manufacture date, expiration date and batch number of the product: SEE LABEL.

CONTRAINDICATIONS

ADVERSE FEFECTS

a) Metal sensitivity Reactions in patients have rarely been informed. Implantation of strange material in tissues results in histological reactions involving various sizes of macrophages and fibroblasts. The importance of this chemical effect is unclear, as well as similar changes may occur as a precursor, or during the healing process; b)- Delay or no consolidation of bone that can lead to breakage of the implant; c) - Mechanical loosening, twisting, dismemberment, breach of the components - may be the result of defective fixing or occult infection; d) - superficial and/or profound infection; e)- Shortening of the segment and the spinal column due to bone resorption; f)- Pain, discomfort or abnormal sensations due to metal implant component placement or due to migration of implantable mechanical components, and thrombosis; g) - Pseudo-arthrosis; h)-Fracture; i)- Paralysis; j)- Hematomas and bleeding in blood vessels; k)- Pulmonary embolism; m) - Neural or neurological problems caused by surgical trauma such as: gastrointestinal disorders; n) - Early or late failure of the implantable meta components; o) Chronic pressure sores; p)- Superficial or muscular sensitivity in patients with inadequate tissue coverage in the operated area; q)- Difficulty in performing physical activities; r)- Allergy; s)- Illness and Death: in all surgical procedures there is an incidence of disease and death. The patient must be warned (by the surgeon) about this possible incidence before undergoing the surgery.

PATIENT INFORMATION

THE INFORMATION ABOUT THE FOLLOWING ITEMS MUST BE SUPPLIED TO PATIENTS:

Indications - Contraindications - Information Use Policies - Possible Adverse Effects - Warnings and Precautions.

The patient must be informed about the importance of post surgery monitoring. The lack of monitoring spoils the detection of post-surgical problems such as loosening of components or occurrence of osteolysis. The non-realization of revision surgery when the loosening of components or osteolysis occurs may result in progressive loss of periprosthetic bone stock.

The patient must be informed about the type of implant used, its brand, code, name, batch, the surgeon's name and surgery

Performance foreseen in General Requirements of ANVISA Regulation, which provides the essential safety and efficacy requirements for products, described in the resolution - RDC 56, April 6th 2001:

The non-expansive Cages - BM are manufactured in Titanium Alloy, according to ASTM F136 proven by laboratory testing reports and quality certification of the product. These materials are acceptable for the manufacture of implants as defined in Annex A of the NBRISO21534 standard. The manufacture and packaging of non-expansive Cages – BM are based on

The non-expansive Cages – BM are distributed on a non-sterile form and packed in transparent and nontoxic PVC plastic envelope, sealed in machines assuring its integrity. The non-expansive Cages – BM are distributed with sufficient information focusing on indications and sufficient information for the product correct handling as well as the safety of the ones involved with its handling.

ALLOWABLE COMBINATIONS WITH OTHER MATERIALS

Metallic alloy Implants of different chemical compositions are from different manufacturers are not recommended for reasons of chemical, physical, biological and functional incompatibility.

The following related ancillary components must be purchased separately, because they are not parts of this product. - Spine Fixation System - Pedicular - ANVISA Registration Number 80128580081 (neither an object of this register, nor part of Implant-System for Posterior Spine Setting - BM Plate - ANVISA Registration Number 80128580083 (neither an object of

this register, nor part of this product)

Bone Graft (neither an object of this register nor part of this product): The implants may be filled with any bone graft. Such choice is the surgeon's responsibility. This bone graft begins to grow through the existing holes in the walls of the cages forming a solid bone (fusion) which shall hold the vertebrae together. This process is known as intervertebral

ACCESSORIES

e non-expansive Cages – BM do not have any accessory with the aim of integrating the medical product.

IMPLANTED PRODUCT EVALUATIONS

After the implantation during the surgery, the concerning professional must execute radiological control to verify the correct positioning of the product. The professional must be in charge of leading clinical and radiological evaluations after the surgery, since it is his/her own responsibility to do so. It has to be done according to his/her own decision toward frequency in order to check the implant condition and bone consolidation evolution. In case the product is found out of the correct positioning, or indicate any non-compliance, it is the surgeon responsibility to take the most appropriate corrective action.

USEFUL INFORMATION FOR PREVENTING RISKS DERIVING FROM IMPLANTATION

ation contraindications, instructions for use and all information contained in the "Instructions for To reduce the risks of implant Use" of the product must be strictly followed.

PRODUCT DISPOSAL

The non-expansive Cages – BM explanted from patients must be properly disposed by the hospital. It is the hospital responsibility the complete adulteration of the implant avoiding its reuse. Biomecanica recommends that the explanted implants are deformed mechanically by using a hammer or an impact press. Then, it must be properly identified with the As described in Resolution No. 2605 of August8th 2006, implantable devices of any kind classified as single-use, are

ADDITIONAL INFORMATION ABOUT THE PROCEDURE BEFORE USING THE PRODUCT, STERILIZATION METHOD

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the non-expansive Cages – BM metal components get in contact with tissue and body fluids. Explanted products must be treated as

DECONTAMINATION, CLEANING AND STERILIZATION The product comes previously cleaned in its original packaging, and ready for sterilization. It must be handled with care to avoid contamination. If the product is subjected to conditions that undermine their state of cleanliness, the implant must be cleaned before sterilization. Below we suggest some procedures for cleaning:

The cleaning can be developed through manual or mechanical methods

EQUIPMENT FOR CLEANING:

ultrasonic Washe
 sterilizer Washer

Manual Cleaning: the procedure is performed manually, where the dirt is removed through physical action with aid of detergent, water and artifacts such as a sponge and a brush. The manual cleaning is strongly recommended because it causes less damage to surgical implants and instruments. In manual cleaning, water must be used at ambient temperature, nylon brush should be used instead of steel or abrasive one because it can damage the material protective layer. Decontamination should be done by using a bactericide and broad spectrum antifungal solution, avoiding aggressive cleaning products, brushes made of metal so that the quality of metallic implants are not damaged. In case of using cleaning machines, the implants and surgical instruments must be placed in drawers so that they do not bump each other avoiding damage to the material.

each other avoiding damage to the material.

The ones in charge of cleaning the implants and surgical instruments must be careful about the type of cleaning product used respecting the concentrations specified by the manufacturer. The implants and surgical instruments must be thoroughly rinsed until all residues are removed and they must be dried immediately after cleaning.

TYPES OF CLEANERS

Enzymatic cleaners: they are basically composed by enzymes, surfactants and solubilising. A balanced combination of these elements means that product can remove the organic matter of the material in a short period of time.

Enzymatic solutions: have excellent action for cleaning, but do not have bactericidal and bacteriostatic activity.

Enzymes: are substances produced by living cells and chemical reactions that govern the process. Once produced by cells,

an enzyme can be isolated and will maintain its catalytic properties, if certain conditions are maintained in its manufacture Enzymes are classified into three major functional groups depending on the type of substrate that they will affect: protease, lipase and amylase that act on protein substrates, fats and carbohydrates, which tend to solubilize themselves and come off the items. Currently, cleaning of articles of complex configuration is recommended to ensure cleanliness.

DISINFECTION:

It is a process that destroys micro-organisms, pathogenic or not, of the items, except for bacterial spores, by physical or

LEVELS OF DISINFECTION:

High level: destroys all microorganisms except the high number of spores -2% Glutaraldehyde - 20-30 minutes.

Indication: preferably in hospitals.

Middle level: eliminates vegetative bacteria, most viruses, fungi and mycrobacteria

-1% sodium hypochlorite -30 minutes.
Indication: For UBS, day care, rest homes, nursing homes.
Low level: eliminates most bacteria, fungi and some viruses, but does not eliminate mycrobacteria

- Sodium hypochlorite 0,025%

This product supplied is not sterile. Before using, it should be sterilized. We recommend the sterilization steam autoclave in the hospital (ISO 11134: 1994 Sterilization of health care products - Requirements for validation and routine control

It is an apparatus (as for sterilizing) using superheated steam under high pressure

OPERATING INSTRUCTIONS

To start the equipment, check: if the circuit breaker is connected if the waterworks control device is open;

1 - open the equipment door: 2 - accommodate the material to be sterilized properly:

5 - connect the electricity key; 6 - the cycle will run automatically in sequence

7 – when the light goes on it indicates "end of cycle", partially open the door for ten (10) minutes, in order to cool the material.

It is recommended that they apply the following parameters of physical sterilization in an autoclave (saturated steam):

Time of Exposure Temperature

Conventional (1 atm pression)	121°C (250°F)	30 minutes
Conventional (1 atm pression)	132°C (270°F)	15 minutes
Gravity	132°C (270°F)	45 minutes
Vacuous	132°C (270°F)	7 minutes

Note: The time should be marked when the heat of the sterilization chamber reaches the desired temperature.

FOR BEST CLARIFICATION READ THE INSTRUCTION MANUAL OF THE AUTOCLAVE Another method of sterilization: the following sterilization method can be used instead of the autoclave as defined by the

Sterilization by ethylene oxide (ETO) - parameters and procedures established in the protocol for validation and EN550;

1994 Sterilization of Medical Devices - Validation and Routine Control of Ethylene Oxide Ste

CAGE NO EXPANSIVO

CARE WITH STERILIZED PRODUCTS

CARE WITH STERILIZED PRODUCTS

Conditions for storage of sterile items related to the environment; it must be clean, airy and dry; it should be restricted to staff; related to the product: after the sterilization process do not place it on cold surface (stone or stainless steel), use baskets or containers with small openings until it cools down; Mantle (unbleached cotton fabric, nonwoven, surgical grade paper, creped paper, film paper, tyve or perforated metal boxes) should remain intact and not be handled to prevent packages from tearing or releasing the seal; they should be stored in locked cabinets with shelves; the shelves must be identified in order to facilitate the removal of the material; the material should be stored according to the expiration date of the sterilization to facilitate distribution and not to keep products that are not valid in stock; they should be stored separately from non-sterile products to

MAILING OF MATERIAL TO BE ANALYZE BY THE MANUFACTURER

The implants to be analyzed by the manufacturer must be cleaned in the hospital using a bactericide and broad spectrum antifungal solution. Then, they should be disinfected or sterilized in a steam auto clave or using ethylene oxide. They should be sent to Biomecânica in a packaging identified with the cleaning and sterilization method, and data of the product.

AFTER SALE SUPPORT (CLIENT COMPLAINT)

If there is a need to make a complaint of non expansive cages - BM related to any adverse effect that affects the safety of the user, such as if the product is not working, implantable metal component damage, death or serious problems related to these components, the responsible surgeon must communicate this adverse event for the national health agency and Biomecanica through e-mail or phone sac@biomecanica.com.br +55 14 2104 7926. In cases of doubt, the responsible surgeon or the health professional can make the communication of adverse events through the notification system in the Health Surveillance

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Sistema da Qualidade Certificado

Quality System Certified Quality System Certified ISO 9001 ISO 13485 REV01

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